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**A GUIDE TO THE
GENERAL
PROPERTIES OF
THE RESIN USED IN
KARTELL LABWARE**

The resins can be split up into two main groups:

Thermosetting resins

Are formed by heat. During moulding operations they complete their space-lattice achieving a rigid irreversible structure.

Thermoplastic resins

Are meltable resins with linear structure instead of lattice. During moulding operations they are transformed without a chemical modification into objects thus giving a reversible formation process, that can be repeated without changing their thermoplastic properties. Thermoplastic resins are the materials commonly used in Kartell plastic labware production. For this reason we give below a general description of these resins underlining their structural, mechanical, chemical and physical properties.

Polyolefins

Polyolefins are high molecular weight hydrocarbons. They include low-density, linear low-density and high-density polyethylene, polyallomer, polypropylene and polymethylpentene. All are break-resistant, non-toxic, and non-contaminating. These are the only plastics lighter than water. They easily withstand exposure to nearly all chemicals at room temperature for up to 24 hours. Strong oxidizing agents eventually cause embrittlement. All polyolefins can be damaged by long exposure to ultraviolet light.

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Polyethylene (PE)

The polymerisation of ethylene results in an essentially straight-chained, high molecular weight hydrocarbon. Branching (side chain formation) occurs to some extent and can be controlled.

Like other polyolefins, the polyethylenes are chemically inert. Strong oxidizing agents will eventually cause oxidation and embrittlement. They have no known solvent at room temperature. Aggressive solvents will cause softening or swelling, but these effects are normally reversible.

Low-density polyethylene (LDPE) has more extensive branching, resulting in a less compact molecular structure. High-density polyethylene (HDPE) has minimal branching, which makes it more rigid and less permeable than LDPE. Linear low-density polyethylene (LLDPE) combines the toughness of low-density polyethylene with the rigidity of high-density polyethylene.

Polypropylene (PP)

It is similar to polyethylene, but each unit of the chain has a methyl group attached. It is translucent, autoclavable, and has no known solvent at room temperature. It is slightly more susceptible to strong oxidizing agents than conventional polyethylene because of its many branches (methyl groups, in this case).

Polymethylpentene (PMP or TPX)

It is similar to polypropylene, but it has an isobutyl group instead of a methyl group attached to each monomer group of the chain. Its chemical resistance is close to that of PP. It is more easily softened by some hydrocarbons and chlorinated solvents. PMP is slightly more susceptible than PP to attack by oxidizing agents. Its excellent transparency, rigidity, and resistance to chemicals and high temperatures make PMP a superior material for labware. PMP withstands repeated autoclaving, even at 150°C. It can withstand intermittent exposure to temperatures as high as 175°C. Products made of polymethylpentene are brittle at ambient temperature and may crack or break if dropped from benchtop height.

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


Polystyrene (PS)

Rigid and non-toxic, with excellent dimensional stability and good chemical resistance to aqueous solution but limited resistance to solvents. This glass-clear material is commonly used for disposable laboratory products. Products made of polystyrene are brittle at ambient temperature and may crack or break if dropped from benchtop height.

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Polyvinyl Chloride (PVC)

It is similar in structure to polyethylene but each unit contains a chlorine atom. The chlorine atom renders it vulnerable to some solvents, but also makes it more resistant in many applications (PVC has extremely good resistance to oils and very low permeability to most gases). Polyvinyl chloride is transparent and has a slight bluish cast. When blended with phthalate esters plasticizers, PVC becomes soft and pliable, providing tubings of any dimension.

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Acetals

Due to their chemical composition, regular and highly crystalline structure, these resins have physical properties that cannot be obtained either with metals or with other plastic material. In addition these resins have high mechanical resistance and stiffness together with low coefficient of friction and high wear resistance.

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Polyoxymethylene (POM)

It is produced by polymerisation of formaldehyde. Acetal retains its dimensions and other properties at elevated temperatures. It offers excellent resistance to strong acids and bases. Naturally opaque.

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Polyamidic resins


Nylon (PA6) .

This is a group of linear polymers with repeated amide linkages along the backbone. These are produced by an amidation of diamines with dibasic acids, or polymerisation of amino acids. Nylon is strong and tough. It resists abrasion, fatigue and impact. Nylon offers excellent chemical resistance with negligible permeation rates when used with organic solvents. However, it has poor resistance to strong mineral acids, oxidizing agents and certain salts.

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Polycarbonates (PC)

This is a special type of polyester, in which dihydric phenols are joined through carbonate linkages (O-CO-O). These linkages are subject to chemical reaction with bases, concentrated acids, etc., and make PC soluble in various organic solvents. PC is window-clear, amazingly strong and rigid. It is autoclavable, non-toxic, and the toughest of all thermoplastics. PC maintains its resistance to impact in a wide range of temperatures and even under very severe environmental conditions. It withstands both low and high temperature from -50°C up to +130°C and has extremely good optical properties together with a high resistance to sunlight exposure (UV radiation).

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
Acrylic resins (PMMA)

Rigid, transparent, resistant to the atmospheric agents, it can replace glass in all the applications where temperature is below 90°-100°C.

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Fluorocarbons

Typical fluorocarbons are Polytetrafluoroethylene (PTFE), Fluorinated Ethylene Propylene (FEP), perfluoroalkoxy (PFA), Polyvinylidene Fluoride (PVDF), Ethylene-tetrafluoroethylene (ETFE) and Ethylene-Chlorotrifluoroethylene (E-CTFE).

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

Polytetrafluoroethylene (PTFE)

PTFE is a high molecular weight polymer with the following fundamental characteristics: chemical resistance to all reagents and solvents, exceptional heat stability (nonflammable), superior self-lubrication and high impact characteristics, toughness and flexibility even at low temperatures. Since 1940 PTFE's have been extensively used in industry for gaskets, bearings, «O» rings etc. They can withstand temperatures from -270°C to +250°C (up to 300°C for short periods).

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Tetrafluoroethylene -

perfluoropropylene (FEP)


FEP is translucent, flexible and feels heavy because of its high density. It resists all known chemicals except molten alkali metals, elemental fluorine and fluorine precursors at elevated temperatures. It should not be used with concentrated perchloric acid. FEP withstands temperatures from -200°C to $+205^{\circ}\text{C}$, and may be sterilised repeatedly by all known chemical and thermal methods. It can even be boiled in nitric acid.

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
Perfluoroalkoxy (PFA)

PFA is translucent and slightly flexible. It has the widest temperature range of the fluoropolymers - from -270°C to $+260^{\circ}\text{C}$ - with superior chemical resistance across the entire range. Compared to PTFE (at $+277^{\circ}\text{C}$), it has better strength, stiffness and creep resistance. PFA also has a low coefficient of friction, outstanding non-stick properties and is flame-resistant.

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Polyvinylidene Fluoride (PVDF)

This is a fluoropolymer with alternating CH_2 and CF_2 groups. PVDF is an opaque white resin. Extremely pure, it is superior for non-contaminating applications. Mechanical strength and abrasion resistance are high, similar to ECTFE. It resists UV radiation. The maximum service temperature for rotationally-moulded PVDF tanks is 100°C . Up to this temperature, PVDF has excellent chemical resistance to weak bases and salts, strong acids, liquid halogens, strong oxidizing agents and aromatic, halogenated and aliphatic solvents. However, organic bases and short chain ketones, esters and oxygenated solvents will severely attack PVDF at room temperature. Fuming nitric acid and concentrated sulphuric acid will cause softening.

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Ethylene - tetrafluoroethylene (ETFE)

This is white, translucent and slightly flexible. It is a close analog of Teflon fluorocarbons, an ethylene tetrafluoroethylene copolymer. ETFE shares the remarkable chemical and temperature resistance of PTFE and FEP, and has even greater mechanical strength and impact resistance.

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**Ethylene -
chlorotrifluoroethylene
(E-CTFE)**

This is an alternating copolymer of ethylene and chlorotrifluoroethylene. This fluoropolymer withstands continuous exposure to extreme temperatures and maintains excellent mechanical properties across this entire range (from cryogenic temperatures to 180°C). It has excellent electrical properties and chemical resistance, having no known solvent at 121°C. It is also non-burning and radiation-resistant. Its ease of processing lends itself to wide range of products.

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**Biological properties of
plastics**

In general, we consider most of the plastics used in Kartell Lab-ware to be biologically inert. For example LDPE, HDPE, PP, PMP (TPX), PC, PS, PTFE, FEP, and PFA have been shown to be non-toxic to cell cultures. Distilled water for preparing culture media is often collected and stored in polyethylene containers.

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**ABBREVIATIONS,
TEMPERATURE AND
CHEMICAL RESISTANCE
OF PLASTICS**

The table below lists commonly used abbreviations for plastics. This list covers plastics commonly employed in the manufacture of plastic laboratory ware.

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Temperatures appearing in parentheses: limits tolerated for intervals only.

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**LIST OF PLASTICS AND
THEIR CHEMICAL
RESISTANCES TO
SUBSTANCE GROUPS**

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● Good resistance; no, or only minor, damage resulting from exposures of more than 30 days

● Marginal resistance; for some types of plastics, extended exposure can result in damage (hairline cracks, loss of mechanical strength, discolouration, etc.)

● Non resistant; exposure can lead to deformation or destruction.
- Résistance excellente.

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Chemicals est
20° or 50°C, resp.

		ETFE E-CTFE 50	PA6 20	PC 20	PE 50	PMMA 50	POM 50	PP 50	PS 50	FEP PFA PTFE 50	PVC 20	SAN 50	PMP (TPX) 50
Acetaldehyde	CH ₃ CHO	●	●	●	●	●	●	●	●	●	●	●	●
Acetone	CH ₃ COCH ₃	●	●	●	●	●	●	●	●	●	●	●	●
Acetophenone	C ₆ H ₅ COCH ₃	●	●	●	●	●	●	●	●	●	●	●	●
Ethyl acetate	C ₂ H ₅ COOCH ₃	●	●	●	●	●	●	●	●	●	●	●	●
Ethyl alcohol (96%)	C ₂ H ₅ OH	●	●	●	●	●	●	●	●	●	●	●	●
Ethylene chloride	HBC ₂ CHCl ₂	●	●	●	●	●	●	●	●	●	●	●	●
Ethylene glycol	HOH ₂ C-CH ₂ OH	●	●	●	●	●	●	●	●	●	●	●	●
Allyl alcohol	H ₂ C=CH-CH ₂ OH	●	●	●	●	●	●	●	●	●	●	●	●
Aluminium chloride	AlCl ₃	●	●	●	●	●	●	●	●	●	●	●	●
Formic acid 85%	HCO-OH	●	●	●	●	●	●	●	●	●	●	●	●
Ammonia 25% ac. soln		●	●	●	●	●	●	●	●	●	●	●	●
Ammonia	NH ₃	●	●	●	●	●	●	●	●	●	●	●	●
Ammonium chloride ac. soln	NH ₄ Cl	●	●	●	●	●	●	●	●	●	●	●	●
Amyl acetate		●	●	●	●	●	●	●	●	●	●	●	●
(Pentyl acetate)	CH ₃ -COOC ₅ -H ₁₁	●	●	●	●	●	●	●	●	●	●	●	●
Amyl alcohol	C ₅ H ₁₁ OH	●	●	●	●	●	●	●	●	●	●	●	●
Aniline	C ₆ H ₅ NH ₂	●	●	●	●	●	●	●	●	●	●	●	●
Arsenic acid	H ₃ AsO ₄	●	●	●	●	●	●	●	●	●	●	●	●
Benzaldehyde	C ₆ H ₅ -CHO	●	●	●	●	●	●	●	●	●	●	●	●
Benzene		●	●	●	●	●	●	●	●	●	●	●	●
Benzole	C ₆ H ₆	●	●	●	●	●	●	●	●	●	●	●	●
Lead acetate	(CH ₃ -COO) ₂ Pb 3H ₂ O	●	●	●	●	●	●	●	●	●	●	●	●
Lead acetate aq. soln		●	●	●	●	●	●	●	●	●	●	●	●
Boric acid 10% ac soln	H ₃ BO ₃	●	●	●	●	●	●	●	●	●	●	●	●
Hydrobromic acid 69%	HBr-2H ₂ O	●	●	●	●	●	●	●	●	●	●	●	●
Butyl acetate	CH ₃ COOC ₄ H ₉	●	●	●	●	●	●	●	●	●	●	●	●
Calcium chloride aq. soln	CaCl ₂	●	●	●	●	●	●	●	●	●	●	●	●
Calcium hypochloride aq. soln	Ca(OCl) ₂	●	●	●	●	●	●	●	●	●	●	●	●
Chlorine	Cl ₂	●	●	●	●	●	●	●	●	●	●	●	●
Chlorobenzene	C ₆ H ₅ Cl	●	●	●	●	●	●	●	●	●	●	●	●
Chloroform	CHCl ₃	●	●	●	●	●	●	●	●	●	●	●	●
Chloruft (Chlorine-) water	(HCl)-(HOCl)	●	●	●	●	●	●	●	●	●	●	●	●
Chromic acid 20%	CrO ₃	●	●	●	●	●	●	●	●	●	●	●	●
Chromatosulphuric acid conc.	H ₂ SO ₄ -K ₂ Cr ₂ O ₇	●	●	●	●	●	●	●	●	●	●	●	●
Decahydronaphthalene (Decalini)	C ₁₀ H ₁₈	●	●	●	●	●	●	●	●	●	●	●	●
(Perhydronaphthalene)		●	●	●	●	●	●	●	●	●	●	●	●
Diethyl ether	H ₅ C ₂ -O-C ₂ H ₅	●	●	●	●	●	●	●	●	●	●	●	●
Dibutyl phthalate	C ₈ H ₁₈ (COOC ₄ H ₉) ₂	●	●	●	●	●	●	●	●	●	●	●	●
1.4 Dioxane	C ₄ H ₈ O ₂	●	●	●	●	●	●	●	●	●	●	●	●
Chrystalline (glacial acetic acid)	CH ₃ COOH	●	●	●	●	●	●	●	●	●	●	●	●
Fluorine	F(FI)	●	●	●	●	●	●	●	●	●	●	●	●
Fluorinated hydrocarbons		●	●	●	●	●	●	●	●	●	●	●	●
Hydrofluoric acid 35%	(HF) ₂	●	●	●	●	●	●	●	●	●	●	●	●
Formaldehyde	HCHO	●	●	●	●	●	●	●	●	●	●	●	●
Glycerine	CH ₂ OH-CHOH-CH ₂ OH	●	●	●	●	●	●	●	●	●	●	●	●
Erea(carbamide and other)	H ₂ N-CO-NH ₂	●	●	●	●	●	●	●	●	●	●	●	●
Hexane	CH ₃ -(CH ₂) ₄ -CH ₃	●	●	●	●	●	●	●	●	●	●	●	●
Iodine (tincture of.)		●	●	●	●	●	●	●	●	●	●	●	●
Potassium chloride aq. soln	KCl	●	●	●	●	●	●	●	●	●	●	●	●
Potassium hydroxide	KOH	●	●	●	●	●	●	●	●	●	●	●	●
Potassium permanganete aq. soln	KMnO ₄	●	●	●	●	●	●	●	●	●	●	●	●
Aqua regia	(HNO ₃) (HCl)	●	●	●	●	●	●	●	●	●	●	●	●
	HNO ₃ +3HCl NOCl + Cl ₂ + 2H ₂ O	●	●	●	●	●	●	●	●	●	●	●	●
Cupric sulphate aq. soln	CuSO ₄	●	●	●	●	●	●	●	●	●	●	●	●
Magnesium chloride aq. soln	MgCl ₂	●	●	●	●	●	●	●	●	●	●	●	●
Methyl alcohol (methanol)	CH ₃ OH	●	●	●	●	●	●	●	●	●	●	●	●
Methylene chloride	CH ₂ Cl ₂	●	●	●	●	●	●	●	●	●	●	●	●
Sodium carbonate	Na ₂ CO ₃	●	●	●	●	●	●	●	●	●	●	●	●
Sodium dichromate	Na ₂ Cr ₂ O ₇	●	●	●	●	●	●	●	●	●	●	●	●
Sodium hydroxide	NaOH	●	●	●	●	●	●	●	●	●	●	●	●
Oxalic acid, 10% aq. soln	HOOC-COOH	●	●	●	●	●	●	●	●	●	●	●	●
Ozone <0,5 ppm (active oxygen)	O ₃	●	●	●	●	●	●	●	●	●	●	●	●
Perchloroethylene	Cl ₂ C=CCl ₂	●	●	●	●	●	●	●	●	●	●	●	●
Phenol 100% (crystals)	CH ₃ H ₃ OH	●	●	●	●	●	●	●	●	●	●	●	●
Phosphoric acid	H ₃ PO ₄	●	●	●	●	●	●	●	●	●	●	●	●
Phosphorus trichloride	PCl ₃	●	●	●	●	●	●	●	●	●	●	●	●
Pyridine	C ₅ H ₅ N	●	●	●	●	●	●	●	●	●	●	●	●
Mercury	Hg	●	●	●	●	●	●	●	●	●	●	●	●
Mercurous + mercuric chlorid	Hg ₂ Cl ₂ +HgCl ₂	●	●	●	●	●	●	●	●	●	●	●	●
Nitric acid 50% conc.	HNO ₃	●	●	●	●	●	●	●	●	●	●	●	●
Hydrochloric acid		●	●	●	●	●	●	●	●	●	●	●	●
(hydrogen chloride)	HCl	●	●	●	●	●	●	●	●	●	●	●	●
Sulphuric acid 95%	H ₂ SO ₄	●	●	●	●	●	●	●	●	●	●	●	●
Silver nitrate	AgNO ₃	●	●	●	●	●	●	●	●	●	●	●	●
Carbon tetrachloride	CCl ₄	●	●	●	●	●	●	●	●	●	●	●	●
Tetrahydrofuran	THF	●	●	●	●	●	●	●	●	●	●	●	●
Toluene	C ₆ H ₅ CH ₃	●	●	●	●	●	●	●	●	●	●	●	●
Trichloroethylene	C ₂ H ₃ Cl ₃	●	●	●	●	●	●	●	●	●	●	●	●
Tri-sodium phosphate	Na ₃ PO ₄	●	●	●	●	●	●	●	●	●	●	●	●
Hydrogen peroxide 30%	H ₂ O ₂	●	●	●	●	●	●	●	●	●	●	●	●
Xylene	C ₆ H ₄ (CH ₃) ₂	●	●	●	●	●	●	●	●	●	●	●	●
Zinc chloride 10% aq. soln	ZnCl ₂	●	●	●	●	●	●	●	●	●	●	●	●
Zinc sulphate 10% aq. soln	ZnSO ₄	●	●	●	●	●	●	●	●	●	●	●	●

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STERILISING PLASTIC LABORATORY WARE

Before sterilising any items of plastic laboratory ware, verify that no contamination or residues are present. Their presence could destroy plastics during sterilisation or autoclaving.

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*Autoclavable at 121°C for a period of 20 minutes.

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Observe the tolerated temperature ranges of plastic when autoclaving plastic laboratory ware. Remove any stoppers, fittings, or caps from plastic laboratory ware prior to autoclaving. Plastic vessels should be autoclaved separately from their closures and other fitting. Autoclaving vessels with their closures in place can lead to deformation and destruction of the vessels.

All statements are advisory only, and imply no liability on the part of Kartell. All statements relating to the resistances of plastic laboratory ware to high temperatures, chemicals, and to sterilisation and cleaning procedures have been cautiously formulated, based on statements of raw material manufacturers, on statements appearing in the literature, plus experience gained in actual practice.

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CLEANING PLASTIC LABORATORY WARE

All polyolefins, such as LDPE, HDPE, PP and PMP (TPX), as well as the fluorinated hydrocarbons PTFE, PFA, FEP, ETFE and E-CTFE, have wettable surfaces that are both highly resistant to high temperatures and chemical attack and easy to clean. Slight contamination can be removed using a chemically neutral (pH 7) cleaning agent. Heavy contamination can be removed using an alkaline (pH up to 12) cleaning agent.

We recommend the use of Kartell Cleanilab LM1, LM2 and LA2 cleaning agents. Use only chemically neutral (pH7) cleaning agents on polycarbonate (PC) or polystyrene (PS).

Laboratory ware used in trace analyses should be cleaned in a 1-N hydrochloric acid (HCl) solution for periods of not more than 6 hours and then rinsed in distilled water in order to preclude contamination by cations or anions.

Never use scouring powders or abrasive sponges when cleaning plastic laboratory ware. Use no alkaline cleaning agents on polycarbonate (PC) laboratory ware.

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