## Fax to +49 (0)3328/437-444

## Ultrasonic Device UP400S (400 watts, 24hHz) www.hielscher.com/up400s

Stand- mounted Ultrasonic Device UP400S		
	UIP400S for beaker sonication in sound protection box	UP400S with flow cell for continuous ultrasonication
Technical Details	UP400S: Powerful ultrasonic lab device (400W, 24kHz) for stand- mounted operation, automatic frequency tuning, amplitude adjustable from 20 to 100%, pulse adjustable from 0 to 100%, dry running protected, with 9-pin DSUB interface, in portable case, IP40 grade, titanium horn	
Applications	Homogenization, Dispersing & Deagglomeration (e.g. nano particles), Wet-milling & Grinding (e.g. nano materials), Emulsifying, Lysis & Cell Disintegration, Extraction, Sonochemistry	
Advantages	<ul> <li>Effectiveness: High-power lab device for sonication of larger volumes from 5 mL up to 4L in the batch and approx. 10 to 50 litres per hour with flow cell.</li> <li>Flexibility: Beaker sonication and continuous sonication for applications in lab &amp; bench-top scale and for small production; various accessories allow for the adaption to specific applications.</li> <li>Simplicity: The UP400S is simple to use and easy to clean.</li> <li>Accessories such as a timer allows for time controlled operation and a PC-control make power and temperature controlling easy.</li> </ul>	
Items: □ Ultrasonic Processor UP400S ∘ 115 volt ∘ 230 volt □ Sonotrode H7 with 7mm tip diameter □ Sonotrode H22 with 22mm tip diameter □ Stand ST1-16 □ ST1-Clamp □ Sound Protection Box SB1-16 □ Power Meter ∘ 115 volt ∘ 230 volt		
ielscher ultrasonics		

Hielscher Ultrasonics GmbH www.hielscher.com; info@hielscher.com Tel.: +49 (0)3328/437-3 Fax: +49 (0)3328/437-444

MANEKO s.r.o., Na Pískách 71, 160 00 Praha 6; tel.: 233 335 638 - 9, fax: 233 332 656; maneko@maneko.cz; www.maneko.cz