

Pure silence, advanced ergonomics **BIOWIZARD** GoldenLINE





The superior advantages of

BIOWIZARD GoldenLine

Microbiological Safety Cabinet class || certified according to EN 12469





Silent working conditions as a standard feature

Save your brain, concentrate on work

European Agency for Safety and Health at work has made work related noise one of its main themes. Current Microbiological Safety Cabinets have high noise levels and the EN norm settles for under 65 dB(A).

Noise is undesired and unpleasant for workers. It puts load on brains and can harm physically and mentally.

High sound level not only affects hearing, stress level and blood pressure but also weakens concentration, creativity and comprehension.

Kojair Tech Oy has developed a Microbiological Safety Cabinet according to EN12469 with sound levels from 48 dB(A).

Practical ergonomics

- Extremely silent, comfortable operation
- The electric sliding front panel closes completely, also available in a manual version
- Efficient light tubes ensure shadowless lighting over the whole working area. Max. >2000 lux.
- Laminated safety glass side windows
- Comfortable arm support made of material that feels warm against the arm
- The backwards-slanted front side gives a clear view through the sliding door and guarantees a comfortable working position

Silence

- Most silent on the market
- Under 52 db(A) *

Lights

- Adjustable lights
- Over 2000 lux *
- No shadows
- Glare free

* in Biowizard Golden 130

Ergonomics

- Comfortable working position
- Easy-to-use control panel
- Easy to clean
- Adjustable work height (700 - 900 mm) can be chosen on site

Service and maintenance

• All maintenance work via front panel

Design

- Nice lines, good looking curves
- Fresh colour

Equipment

Comes fully loaded as standard



Adjustable lights

- Light exceeds EN 12469
- Lights fully adjustable between 750 and 2000 lux
- No shadows
- Anti-glare shields
- Light tubes can be changed from outside







Easy to clean

- Easily cleaned front glass
- All surfaces are smooth without any screws or crevices
- The polished, standard stainless steel three-piece table top is easy to remove for cleaning
- Stainless steel bottom sink with round corners
- Window opens fully for easy inside cleaning
- The cabinet is easy to clean with normal cleaning solvents

BIOWIZARD

Technical details and options



365

100	130	170	200
975	1280	1625	1930
885	1190	1535	1840
600	600	600	600
720	720	660	660
785	1030	1310	1550
285	380	495	595
1070	1410	1805	2145
200	200	250	250
487	640	812	965
2x55	2x55	2x55	2x55
180	240	290	330
740	940	940	1170
350	400	450	600
	885 600 720 785 285 1070 200 487 2x55 180 3480	975 1280 975 1280 885 1190 600 600 720 720 785 1030 285 380 1070 1410 200 200 487 640 2x55 2x55 180 240	No. No. 975 1280 1625 885 1190 1535 600 600 600 720 720 660 720 720 660 785 1030 1310 285 380 495 1070 1410 1805 200 200 250 487 640 812 2x55 2x55 2x55 180 240 290 740 940 940

*) Possible also for 220-240 V AC, 60 Hz

Transport height with base 1990 mm, without base 1440 mm

Technical details

- Adjustable base stand
- Slanted electrical front window
- Constant measurement of the airflow velocities
- Magnetic safety valve for gas
- Three piece table top
 - Smooth touch user interface with a real time clock
 - 1/2 speed switch
 - Media valve installation channels
- PCB interface for service info
- UV usage hour center
- Socket / UV socket timers

- Options
- Exhaust connection fitting
- Clean room compatible electrical base stand
- One piece table top with lowered middle section
- UV light with magnetic fastening
- Fully integrated double exhaust filtration
- Anti flowback damper
- Potential free connection
- Extra double sockets, compatible with different European standards
- Extra gas and other media valves
- Foot pedal for front window control

Ask for customized applications for your special needs; microscopes, IVF, weighing, isotope solutions etc.

Service

- Access to the follow-up information via the front panel
- Comprehensive Kojair training system
- Spare parts available, guaranteed for 10 years (filters, lights etc.)
- Installation and commissioning
- Easy performance of annual testing procedures, manuals, specs etc.
- All maintenance can be performed through the front
- Due to enlargened HEPA- filter media area, filter life time is extra long

Service measurements with versatile, calibrated instruments









Balometer



10

KI discus tester



Illumination meter

5/2009 Pictorino Oy Kirjapaino Hermes Oy 3.000

Leading expertise in clean air technology



Kojair Tech Oy Teollisuustie 3 FI-35700 Vilppula Finland Tel. +358 3 471 7000 Fax +358 3 471 7077 sales@kojair.com



