### standards and guidelines

Designed and manufactured to comply with the following standards and auidelines:

European Medical Device Directive 93/42/EEC

89/336/EEC European Electromagnetic Compatibility Directive

EN 13060 Small Steam Sterilizers IEC 61010 Equipment Electrical Safety ASME VIII div 1 Pressure Vessel design HTM 2031 Clean Steam for Sterilization

\* With air leak detector test kit fitted (Part code: 309042)

ISO 9001 Quality Management Systems

ISO 13485 Medical Device Quality Management Systems

#### quarantee

Pressure chamber and heating element 10 years









Prestige Medical Limited East House, Duttons Way, Shadsworth Business Park,

Tel: +44 (0) 1254 682 622 Fax: +44 (0) 1254 682 606 www.prestigemedical.co.uk sales@prestigemedical.co.uk

Registered in England Reg No. 2826793 Vat No. 633707835

Registered Office: Unit 1, First Avenue Maybrook Industrial Estate Minworth, Sutton Coldfield West Midlands, B76 1BA

Prestige Medical is a registered trademark



MANEKO spol. s r.o. Na Pískách 71 160 00 Praha 6

Tel.: +420 233 335 638 - 9 Fax: +420 233 332 656

E-mail: maneko@maneko.cz

www.maneko.cz



Due to a continuous process of development, Prestige Medical reserves the right to amend the specification of products and services without notice.

## High capacity for faster throughput





















Take the step to better sterilization and effective infection control



The RCVS Guidelines clearly state "The practice must provide an autoclave, vacuum or non-vacuum or other recognised sterilization system, for the effective sterilization of instruments and equipment."

#### The Anima autoclave sets high standards for speed, efficiency and performance for sterilizing instruments.

Highly efficient, yet compact, its capacity to process 6 full sized trays\* is twice as many as its nearest rivals. So, if you're a busy practice, this capacity to cope with heavy demands saves you valuable time and reduces your operating costs.

\*depending on the type of load

- Vacuum autoclaves are for pouched, hollow, hinged and solid instruments.
- Non-vacuum autoclaves are only suitable for solid instruments. (If you pouch your instruments and put them in a non vacuum autoclave, the machine will not dry them. A vacuum machine enables steam to penetrate all the small gaps of hollow instruments, the hinges of scissors etc. whilst a non vacuum machine cannot do this.)



# Anima - Seamless integration for maximum flexibility



#### More features and benefits...

The Anima autoclave has both a vacuum and a non vacuum cycle on one machine, giving you much more flexibility in use, is available with either 16 or 22 litres capacity chamber and is supplied complete with 3 pouch trays. It also features a single use water system which helps reduce wear and tear on your instruments and a thermal cut out switch which means you don't need to call in an engineer in case of overheating due to accidental overload.

The Anima autoclave is made in the UK by Prestige Medical and is only available direct from the manufacturer. We send our own engineer to install, commission and train you on how to use the machine completely free of charge. Every unit is fully tested before it leaves the factory and is supplied complete with a 12 month warranty\*.

\*Excludes door gasket, internal furniture, air filter and consumables.

**Combines seamlessly** with the Ultra U500H ultrasonic bath ask for details on 01254 685 399

