

Vibratory Sieve Shaker AS 200 control General Information

The analytical sieve shakers of the series AS 200 are used in research & development, quality control of raw materials, interim and finished products as well as in production monitoring. The controllable electromagnetic drive offers an optimal adaption for every product. Sharp fractions are obtained even after very short sieving times.

With its all-digital controls, up to 99 sieving programs and calibration certificate the sieve shaker AS 200 control is indispensable for all users who attach importance to precision and operational convenience and need to comply with the guidelines of the ISO 9001.



Application Examples

cement clinker, chemicals, coffee, construction materials, fertilizers, fillers, flours, grains, metals powders, minerals, nuts, plastics, sand, seeds, soils, washing powder, ...

Product Advantages

- NEW: up to 99 sieving programs
- NEW: low noise operation thanks to optimized control
- NEW: USB connection for use of EasySieve software
- suitable for dry and wet sieving
- · excellent separation efficiency even with short sieving times
- efficient electromagnetic drive
- 3-D throwing motion which ensures optimum use of the open sieve area and lets the sample move equally over the whole sieving surface
- setting of sieve acceleration "g" for comparable and reproducible sieving results worldwide
- sieve stack up to 550 mm height
- free digital adjustment of all process parameters (time, amplitude or sieve acceleration, interval)
- interval operation
- optional sieve software EasySieve for control of the shaker, easy evaluation and documentation of results
- fulfils all criteria for measuring equipment related to ISO 9001
- · calibration certificate
- easy operation, ergonomic design
- maintenance-free



Vibratory Sieve Shaker AS 200 control Features

Applications separation, fractioning, particle size

determination

Field of application agriculture, biology, chemistry /

plastics, construction materials, engineering / electronics, environment / recycling, food, geology / metallurgy, glass /

ceramics, medicine / pharmaceuticals

Feed material powders, bulk materials,

suspensions

Measuring range* $20 \mu m - 25 mm$

Sieving motion throwing motion with angular

momentum

Max. batch / feed capacity 3 kg
Max. number of fractions 11 / 23
Max. mass of sieve stack 6 kg

Amplitude digital, 0.2 - 3 mm

Controlled amplitude yes

Sieve acceleration 1.0 - 15.1 g
Time display digital, 1 - 99 min

Interval operation 1 - 99 s
Storable SOPs 99
Suitable for dry sieving yes
Suitable for wet sieving yes
USB interface yes
Including test certificate / can be yes

calibrated

Suitable sieve diameters 100 mm / 150 mm / 200 mm / 203

mm (8")

Max. height of sieve stack 550 mm

Clamping devices "standard", "comfort", each for wet

and dry sieving

Protection code IP 54 / IP 20

Electrical supply data 100-240 V, 50/60 Hz

Power connection 1-phase

W x H x D 417 x 212 x 384 mm

Net weight ~ 35 kg Standards CE

Please note:

^{*}depending on feed material and instrument configuration/settings



Vibratory Sieve Shaker AS 200 control Videolink



http://www.retsch.com/as200control

Function Principle

All sieve shakers of the series AS 200 work with an electromagnetic drive that is patented by RETSCH (EP 0642844). This drive produces a 3D throwing motion that moves the product to be sieved equally over the whole sieving surface. The advantage: high stress capacity, extremely smooth operation and short sieving times with high separation efficiency.