

InSight

URINALYSIS

LABORATORNÍ POTŘEBY A TECHNIKA

MANEKO

Na Pískách 71, 160 00 Praha 6

Tel. : +420 233 335 638-9

Fax: +420 233 332 656

E-mail: maneko@maneko.cz

www.maneko.cz

Urinalysis is an important tool in disease detection, as well as monitoring and screening animal health. Abnormalities can be indicative of diseases of the urinary system as well as other organ systems, including liver function, acid base status, and carbohydrate metabolism.

MS-11 Urine Strips

Test for:

- Urobilinogen
- Glucose
- Bilirubin
- Ketone
- Specific Gravity
- Leucocytes
- Blood
- pH
- Protein
- Creatinine
- Nitrate

with urine protein/creatinine ratio.

MS-2 Urine Strips

Test for:

- Microalbumin
 - Creatinine
- with urine microalbumin/creatinine ratio.



Urine Protein/Creatinine Ratio (UPC)

Quantitative test for the detection of proteinuria in animals.

Screen patients for the early signs of renal disease and monitor thereafter.



*UPC ratio
result in less
than 2
minutes.*



Use MS strips in conjunction with the MS reader.

Follow us:

Tel: +44 (0)1204 669033

Fax: +44(0)1204 669034

Email: sales@woodleyequipment.com

Web: www.woodleyequipment.com

OLD STATION PARK BUILDINGS, ST. JOHNS STREET, HURWICH, BOLTON, LANCASHIRE, BL6 7NY, UK



WOODLEY
EQUIPMENT COMPANY LTD.